

Vaginal pH Test

CLINICALLY-TESTED VAGINAL INFECTION TEST KIT

Intended Use

Measures vaginal pH, which is important for differentiating between a yeast and bacterial infection.

Test Components

pH test strips

Included in Box

2 individually sealed vaginal pH test strips 2 sterile swabs

2 pH charts to interpret results

Directions

Insert cotton swab into vagina to collect secretions. Apply sample to the pH test pad on the test strip. Compare color of pH test pad to the pH chart (included).



FAST & ACCURATE RESULTS

I Quantity 2ct

I NRC 87651-0760-33 I UPC 787651760339

I Package Dimensions L3.23"x H5.14"x W1.18"

I Weight 0.06 LBS

I Case Pack 24

I Gross Case Weight 1.94 LBS

To place an order for AZO® Products Contact Your Hospital Sales Specialist

ICall 1.800.722.3476

I Email culturellehcp@dsm-firmenich.com

I Visit azoprofessional.com